



*THE ORIGINAL PLATE*

## VibroGym Medical Datasheet

VibroGym Medical is used for treatment and prevention and plays a central role in rehabilitation. Based on the principle of 'whole body vibration' it introduces a new, highly effective training method: VibroTraining®.

VibroGym Medical is produced in the Netherlands to ensure top quality. VibroGym Medical is designed to meet the high demands of therapeutic environments. The complex motor construction consists of two 160W-motors built into the vibration platform. They are steered, controlled and synchronized through electronics by Siemens. Furthermore it is equipped with an additional shielding to avoid interferences with medical devices. The great usability of the VibroGym Medical is ensured by easy to handle and maintenance-free technology, and by easily adjustable training programmes and exercises. VibroGym Medical – a perfect symbiosis of stainless steel, functionality, technology, lasting quality and timeless design.





THE ORIGINAL PLATE

## 2 Years guarantee

For the VibroGym platform we offer a two-year guarantee.

## Coming with the product



### VibroStraps

The upper body straps can be adjusted in length; they stimulate specific groups of muscles.



### VibroMat

A mat is provided for exercises that involve parts of the body, other than the feet, coming into direct contact with the platform. The mat has a damping effect and protects clothing.



### VibroMat Deluxe

The mat is foldable so that users can increase the size of the training area and perform ground exercises at different levels.



### VibroTraining® DVD and CD-ROM

On the enclosed DVD the 36 VibroTraining® exercises are vividly demonstrated. The CD-ROM contains different VibroTraining® Programmes (Introduction, Anti-Aging, Anti-Cellulite, Legs/Bums/Tums, Back Strengthening and Fat Loss) to print by yourself.



### VibroTraining® Poster

The high-quality poster illustrates the 36 VibroTraining® exercises in a self-explanatory way.



### VibroTraining® Academy

In order to be able to optimally use the VibroGym, vibration training is very important. At the primary level, you will learn the basic principles of WBV. You will obtain a VibroInstructor certificate and be able to further specialise at the training institute.

## Accessories



### VibroWheels

Makes the VibroGym suitable for mobile use (storage, group training).




### VibroMasterbox

Central drive system for 1 master VibroGym Professional and up to 8 additional VibroGyms.



### VibroPaymentsystem

Payment system to let the client pay for vibration minutes.



### VibroFloormat

Big floormat 100 x 200 x 30 cm, grey.

## Technical Specifications

Material	stainless steel
Dimensions (l x b x h)	800mm x 800mm x 1400mm
Weight	92kg/197lb
Voltage	115/230 VA
Power consumption	350W
Maximum load	140kg/300lb
Frequency settings	30/35/40/50Hz
Force	2.1 kN/3.9kN
Colour	white RAL 9016 and silver-grey, RAL 9006
Operation	membrane keyboard, start/repeat/stop button
Display	Braille, clear colour coding of start (green) and stop (red) buttons
Certification	CE 3044, Kema as notified body, IFI accredited
Warranty	two-year factory guarantee (conditions apply)